



A048

GYNECOLOGICAL EXAMINATION TABLE

Registration code: SK-CD1



SAIKANG/Medical

Healthier
your life

A048

FEATURES AND HIGHLIGHTS

Technical parameters

▪ External (LxWxH)	1870x600x680 mm
▪ Back-rest Adjustment	0-75°(±10%)

Technical configuration


☒ Drawer	4 pcs
☒ Step Stool	1 pc
☒ Retractable Foot Board	1 pc
☒ Lithotomy Rod	1 set
☒ Stainless Steel Basin	1 pc
☒ Paper Roll Holder	1 pc

RICH PARTS


Big drawers, can accommodate different tools for the examination.




Hidden basin, meet the needs while examination. Can be withdraw while no need.



Paper tool, convenient for the change of paper.



Lithotomy Rod, adjustable in height and different angles.



QUALITY CONTROL

CASE STRUCTURE

Base cabinet structure, thickness 0.8mm ensure the table stable. Load capacity more than 250kg. Easy to clean.



PILLOW

Integrated pillow, comfortable for the patient.



BACK SUPPORT

Mechanic design for back rest lifting, support weight max. 250kg.



LEG SECTION

Retractable auxiliary leg section, meet different height of the patient. Can be withdraw inside.



FOOT STEP

Retractable auxiliary foot step, save space, can store inside while not need.



BACK REST LIFTING

Back left design 75°





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
☎ +86 18021231901 · 🕒 24/7
✉ sales@chinasaikang.com
🌐 www.saikangmedical.com



Company video



General catalogue