




SKH041-3

I.V. STAND



SAIKANG Medical

Healthier
your life

SKH041-3

FEATURES AND HIGHLIGHTS

Technical parameters

- Size 500×490×1310-2170 mm

Structure

Stainless steel telescopic with 4 hooks, plastic knob made of ABS, spring inside the tube. H Shaped feet, made of powder coated square tube.



Hook

Stainless steel two hook, each hook can bear 2kg.



Height Adjustment

Height adjustment by turning the plastic screw.



Base

IV stands can stack together to save space.




Castors

Powder coated feet with 4 Castors, double side, two with brake. diameter: 2 inch.



Color

Color option for feet: white, blue, pink, Tiffany blue, yellow. MOQ 500 pcs for each color.





CONSTANT
INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE
BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
☎ +86 18021231901 · 🕒 24/7
✉ sales@chinasaikang.com
🌐 www.saikangmedical.com



Company video



General catalogue