



SKR-AT588

ANESTHESIA TROLLEY

Registration code: SKM-B



SKR-AT588

FEATURES AND HIGHLIGHTS

Technical parameters

- External size (LxWxH) 588x480x891 mm

Technical configuration

☒ Drawers	3pcs
☒ Casters	4pcs
☒ Trash can	2pcs
☒ Debris basket	1pc

Debris Basket

With 1 stainless steel debris basket on one side.



Caster

4 nylon wheels Ø 100 mm with non-marking tires and plain bearings, with one brakes



Material

Made of the E1 Class HPL board and aluminum alloy frame, resistance of water, acid, alkali, fire and easy clean, not easily polluted by reagents, blood, drugs etc., satisfy the daily disinfectant and infection control with long life.



Guardrail

Aluminum alloy guardrail on three side, anti-corrosion and anti-rust.



Drawer / Worktop

Three drawers made of high quality HPL board, with high strength aluminum handles. Laterally extendable worktop, expansion space can be added for writing



Anesthesia Boxes

With eight pcs anesthesia boxes.



Trash Can

Hidden type trash can shelf, with two trash cans.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
☎ +86 18021231901 · 🕒 24/7
✉ sales@chinasaikang.com
🌐 www.saikangmedical.com