



# SK-EH424

Oxygen  
Concentrator



**SAIKANG** / Medical

LEADING BRAND OF CHINA MEDICAL FURNITURE  
Stock Code: 870098

## TECHNICAL PARAMETERS

- / Oxygen purity: 93±3%
- / Flow range: 0-5L
- / Noise: ≤45dB
- / Input voltage: 220V/100V
- / Output pressure: 30-70kPa
- / Power: 350W
- / Weight: 19kgs
- / Size: 380×320×590mm
- / Max. Flow: 5L/min
- / Continuous working time: 10,000 hours
- / Advanced oil-free compressor
- / with nebulizer function
- / CE, ISO13485 approved by TUV
- / accessory: user manual
- / window spongeHumidifier
- / Oxygen nose tube
- / Oxygen mask (adult)

## PRODUCTS FEATURES

America PSA technology offer nature oxygen  
France Imported molecular sieve bed  
Reliable and durable oil free compressor  
Available 24 hours continuous working  
Self diagnostic system with error code indication


## FUNCTIONS

Power off alarm, Overload protection, High/Low pressure alarm, Temperature alarm, Error code indication, Nebulizer, Oxygen purity alarm.

## ACCESSORIES



Humidification bottle



Nebulizer mask



Nasal oxygen canula



Filter



Sponge filter

## QUALITY CONTROL

### MULTIPLE FLOW OPTIONS

5L, multiple flow options




### OPERATION PANEL

LED operation panel, easy to operate

### USING EFFECT

Low noise, maintenance free




### MOLECULAR SIEVE

Original imported molecular sieve, high efficiency and low energy consumption

### QUALITY

American technology, lasting stability



### COMPRESSOR

48hours working continuously pure copper oil-free compressor, prolong the service life





## INNOVATION OF SAIKANG

LEADING BRAND OF CHINA MEDICAL FURNITURE



Design Innovation



Research Innovation



Manufacturing Innovation



Maintenance Innovation



Service Innovation

### JIANGSU SAIKANG MEDICAL EQUIPMENT CO.,LTD.

24hrs: 0086-18021231901   
[www.saikangmedical.com](http://www.saikangmedical.com)

#### SAIKANG (China Headquarter)

No.35 Lehong Road,Zhangjiagang City,  
Jiangsu Province,China.  
TEL: +86-512-58525850  
FAX: +86-512-58520713  
E-mail: [inquiry@jssaikang.com](mailto:inquiry@jssaikang.com)

#### SAIKANG (USA Branch)

Saikang Medical Technology, Inc  
17800 Castleton St Ste 180,  
City of Industry, Ca 91748  
Tel: +1 949648-2776

#### SAIKANG (UK Branch)

UK Healthward International Limited  
Owner(s) Address-DEPT 111,  
196 High Road Wood Green,  
London, N22 8HH  
United Kingdom  
TEL: +44 2032399738

#### SAIKANG (East Africa Branch)

Saikang Medical Solutions (East Africa) Ltd.  
Viraj Complex Unit 3 Mombasa Road,  
P.O Box 14262-00800, Nairobi.  
TEL: +254 733-704 690

#### SAIKANG (Mid East Branch)

Saikang Medical Supplies (Mid East) LLC  
X22-S03 England Cluster,  
International City, Dubai UAE  
TEL: +20 128 8962497



Catalogue



Website

As the circumstance of printing, the photos might have minor difference on color, tactile from real products.  
As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.